

EN

aurafix
COMPRESSION PRODUCTS




YOUR SUPPORT IN LIFE

2019

COMPRESSION
PRODUCTS

aurafix
COMPRESSION PRODUCTS



C  **MPRESSION
PRODUCTS**



 **MEDICAL COMPRESSION STOCKINGS**



 **LYMPHEDEMA ARMSLEEVES**



 **TRAVEL SOCKS**



 **ANTI EMBOLISM STOCKINGS**





INDEX

■ About Us	5
■ Varicose Vein Formation	7
■ Aurafix Yarn Technology	8
■ Medika Compression Stockings	11
■ Anti Embolism Stockings	17
■ Lymphedema Armsleeves	23
■ Travel Socks	29
■ Silver Socks	35
■ Silver Yarn Technology	43
■ Promotional Items	51





We have been creating value for the medical sector for 16 years with the brands we have established under the roof of TSD Health Products.

We carry on a business in nearly 50 countries and we share our experiences with our patients about orthopedic ready to use products, orthosis, prosthesis, compression and products after liposuction.

Health requires care and attention and does not accept faults. With this awareness, we are performing our business with care and passion, we closely follow the innovations in the sector.

We offer the fastest and most reliable solution to the needs of our patients with our corporate approach that we maintain without compromising our traditions. Today, with all our colleagues under the roof of Aurafix, we are aware of the seriousness of our work and stay the course towards our goals without compromising our principles.

The continuous rise of Aurafix and this sample success are being built in cooperation with our colleagues, patients, dealers, doctors and physiotherapists and we continue to grow.

Türkay Gülhan
CEO




#yoursupportinlife


YOUR SUPPORT IN LIFE

VARICOSE VEIN FORMATION

Deep veins are the main drainage system in the lower extremity venous system. In venous circulation, the veins in our legs have to work in the opposite direction to gravity while sending the blood to the lung. This system is provided by the operation of leg muscles and through the uni-directional valves. However, this mechanism may deteriorate in some individuals over time due to various reasons. This causes the veins to be hyperextended and eventually cause the valves to fail to function. Thus, 'varicose veins' characterized by balloon-like protrusions occur under the skin in the whole leg, especially in the lower part of the leg.

Varicose veins usually pose a problem in terms of aesthetics and appearance, but can also lead to serious health problems when they are not treated. Drugless and painless treatment of varicose veins is the use of compression socks.





Lycra, which has the feature of high elongation, comes in existence from frequent knitting of double-layer polyamide fibers surrounding it. The Lycra, a highly durable, strong and flexible fiber, provides the desired compression, while the polyamide fibers keep moisture transfer and air permeability at the highest level. Our technology that we use ensures that our products are friendly to your skin and long-lasting.





Antibacterial Weaving

Aurafix Compression Products have been produced and tested according to the world standards through developed weaving technology with its anti-bacterial properties for the purpose of eliminating the harmfulness of microorganisms and bacteria.



Non-Perspiring Knitting

The yarns and knitting techniques used in Aurafix Compression Products ensure maximum air permeability and minimizes perspiration by throwing out the moisture.



Double Covered Thread

Thanks to the double-covered yarn produced by the high quality standards of Aurafix Compression Products, the effects that may stem from the yarn remain minimal.



Soft Knitting

The techniques used during the knitting of Aurafix Compression Products allow the products to remain soft and provide a friendly and non-irritating contact to the skin.



Non-Slip Band

Stability of the product is ensured thanks to the applied elastic bands that prevent the drip silicon slip in Aurafix Compression Products and discomfort and uselessness caused by the slip are removed.



Pressure Resistant Heel

The weaving and knitting quality used in Aurafix Compression Products ensures the product life to be longer by showing durability to the intensifying pressure in a maximum level, particularly in the heel area.



YOUR SUPPORT IN LIFE





MEDICAL COMPRESSION STOCKINGS



#yoursupportinlife



■ MEDICAL COMPRESSION STOCKINGS

Aurafix compression stockings with world-class pressure profile and knitting quality are a comfortable and quality solution for the need for compression socks with high quality and durable weaving. Aurafix Compression Socks has been designed in accordance with the standard "TS ENV 12718 Compression Socks". Available in CCL1 (18-21 mmHg), CCL2 (23-32 mmHg), CCL3 (34-46 mmHg), CCL4 (49 mmHg and over) pressures and open-closed toe models. The raw materials used have the certificates of Oeko-Tex 100 and minimize the allergic effects that may arise on the skin.

Compression Stockings apply a reduced compression towards the calf and upper leg, with the external pressure being the highest in the ankle. Thus, they accelerate blood flow, providing a proper pressure distribution on the leg and feet.

INDICATIONS

Light Pressure (18-21 mmHg)

- For prevention of venous diseases,
- For light varicosis not accompanied by edema and light varicosis of pregnancy,
- For diabetes and rheumatic diseases leading to vasodilation and loss of elasticity,
- For treatment after venous surgery.

Medium Pressure (23-32 mmHg)

- In certain varicosis with probability to be edema,
- In post-traumatic swelling,
- During the recovery period of simple ulcers,
- In superficial thrombophlebitis,
- For significant varicosis during pregnancy.

CONTRAINDICATIONS

Absolute Contraindications

- Advanced peripheral arterial occlusive disease.
- Decompensated heart failure.
- Septicophlebitis.
- Phlegmasia coerulea dolens

High Pressure (34-46 mmHg)

- In severe structural and post-thrombotic venous insufficiencies,
- For significant varicosis accompanied by edema,
- In secondary varicosis,
- De'rmatosclerosis,
- In the recovery period of severe and recurrent ulcers (after the ulcer has healed).

Very High Pressure (49 mmHg +)

- For treatment of lymphedema
- For elephantiasis

Relative Contraindications

- Secretory dermatosis.
- Incompatibility with compression stockings material.
- Deterioration of sensitivity in the organs.
- Advanced peripheral neuropathy (Diabetes Mellitus).



Antibacterial Weaving



Non-Perspiring Knitting



Soft Knitting



Double Covered Thread




Non-Slip Band



Pressure Resistant Heel



Anti-bacterial properties of Aurafix Compression Products, produced from high-quality double-covered yarns, have been developed and tested in accordance with the world standards. The yarns and knitting technique ensure maximum air permeability and minimizes perspiration by throwing out the moisture. Thanks to its soft knitting, it provides a friendly and non-irritating contact to the skin. Stability of the product is ensured thanks to the applied elastic bands that prevent the drip silicon slip and discomfort and uselessness caused by the slip are removed. It ensures the product life to be longer by showing durability to the applied pressure particularly from the heel area in a maximum level.







	A-D Knee High Closed Toe	A-D Knee High Open Toe	A-G Thigh High Closed Toe	A-G Thigh High Open Toe	A-T Pantyhose Closed Toe	A-T Pantyhose Open Toe	A-T/U for Pregnant Closed Toe
CCL 1 18-21 mmHg	AV-101 ●	AV-100 ●	AV-103 ●	AV-102 ●	AV-105 ●	AV-104 ●	AV-117 ●
CCL 2 23-32 mmHg	AV-107 ●	AV-106 ●	AV-109 ●	AV-108 ●	AV-111 ●	AV-110 ●	AV-118 ●
CCL 3 34-46 mmHg		AV-112 ●		AV-113 ●		AV-114 ●	

Size	cB	cC	cF	cG
1	18-20 cm	28-33 cm	37-47 cm	43-57 cm
2	20-22 cm	30-37 cm	38-52 cm	46-60 cm
3	22-24 cm	33-40 cm	43-58 cm	49-67 cm
4	24-26 cm	35-43 cm	46-61 cm	52-72 cm
5	26-28 cm	37-45 cm	50-66 cm	55-76 cm
6	28-30 cm	39-46 cm	53-69 cm	60-79 cm
7	30-32 cm	41-48 cm	56-71 cm	64-80 cm

* AD - Knee High
cB

* AG - Thigh High
cB / cC / cF

* AT - Pantyhose & for Pregnant
cB / cC / cG



Caramel

Balck

YOUR SUPPORT IN LIFE





ANTI EMBOLISM STOCKINGS



#yoursupportinlife

ANTI EMBOLISM STOCKINGS

Anti-embolism stockings accelerate blood flow and reduce coagulation by applying constant pressure to the legs. They apply a reduced compression towards the calf and upper leg, with the external pressure being the highest in the ankle. Thus, they increase the blood flow in the deep vein system.

Aurafix Anti-Embolism Stockings have been designed in accordance with the standard "TS ENV 12719 Medical Thrombosis Prophylaxis Socks". There are 3 models as below knee, above knee, above knee belted at the waist. Anti-bacterial weaving in the toe area reduces the risk of fungus. There are anti-slip silicone bands on top of the socks. Thanks to its special weaving having no latex, its compliance with the patient with elastic structure is high.



INDICATIONS

- They decrease the risk of deep vein thrombosis in immobilized patients before and after surgery.
- They increase blood flow rate by providing graded pressure.
- They increase blood flow in other venous disorders such as chronic venous stasis.

CONTRAINDICATIONS

- Regional leg conditions that can be aggravated by product;
 - Cellulite • Dermatitis • Gangrene • Fasciitis, recently skin graphite, • Panniculitis
- Pulmonary edema due to severe arteriosclerosis or other ischemic vascular disease
- Pulmonary edema due to severe edema on the legs or congestive heart disease
- Extreme deformity on the leg



Antibacterial Weaving



Non-Perspiring Knitting



Soft Knitting



Double Covered Thread




Non-Slip Band



Pressure Resistant Heel



Aurafix Anti-Emboli Stockings have been produced from high-quality double-covered yarns, anti-bacterial properties of them have been developed during the weaving and tested in accordance with the world standards. The yarns and knitting technique ensure maximum air permeability and minimizes perspiration by throwing out the moisture. Stability of the product is ensured thanks to the elastic bands that prevent the drip silicon slip and discomfort and uselessness caused by the slip are removed. The inspection (control) hole at the bottom of the sock controls the blood circulation. Pressure resistant and colorful heels provide size separation. Thanks to its soft knitting, it eliminates the risk of pressure sores and ensures the socks life to be longer.





CCL 1
13-18 mmHg


A-D Knee High	A-G Thigh High
AV-99AD	AV-99AG
●	●
	AV-99T
	●




Size	cB	cC	cF
XS	18-20 cm	27-34 cm	37-47 cm
S	20-23 cm	31-37 cm	41-53 cm
M	23-26 cm	35-40 cm	47-58 cm
L	26-29 cm	38-44 cm	52-62 cm
XL	29-32 cm	41-48 cm	58-69 cm

* AD - Knee High
cB / cC

* AG - Thigh High
cB / cC / cF




Color



White

Heel Colors According to Size



Beige Yellow Green Red Blue

XS S M L XL

YOUR SUPPORT IN LIFE





LYPHEDEMA ARMSLEEVES



#yoursupportinlife



LYMPHEDEMA ARMSLEEVES

Lymphedema is a condition which is defined as excessive protein accumulation in interstitial space and has a negative effect on the quality of life of the person and causes restrictions on daily life activities. Lymphedema is classified as primary and secondary lymphedema. Primary lymphedema is associated with the developmental abnormality of lymphatic metastasis (lymph vessel) and lymph nodes. Secondary lymphedema develops due to interruption of lymph as a result of radiotherapy, surgical intervention, trauma, inflammation or neoplastic disorder. Secondary lymphedema is more common than primary lymphedema.

Lymphedema is a chronic condition requiring life-long attention. Control of edema reduces disease progression and soft tissue incidence. Pressotherapy is an important part of the treatment of lymphedema. People with mild lymphedema may reduce bloating by wearing medical pressure armbands and socks. The Aurafix Lymphedema Armsleeves have been designed in accordance with the German standards of "RAL-GZ 387/2 Medical Compression Arm sleeves". 4 different models are available as CCL1 (18-21 mmHg), CCL2 (23-32 mmHg) and CCL3 (34-46 mmHg) with pressure and the ones starting on the wrist and including the thumb.

■ INDICATIONS

- In the treatment of primary or secondary lymphedema,
- In the edema after upper extremity surgery or trauma,
- In the Paget-von-Schrotter syndrome (Lat axilla Phlebothrombosis).

■ CONTRAINDICATIONS

- Vascular occlusion (advanced arteriopathy),
- Uncontrolled advanced heart failure,
- Untreated inflammatory phlebitis,
- Lymphedema that becomes prominent and that does not drop,
- Advanced peripheral neuropathies.

■ PARTIAL CONTRAINDICATIONS

- Immobility (being confined to bed),
- Skin infections, running skin diseases,
- Sensitivity to the material used,
- Sensitivity deterioration in organs.



Antibacterial Weaving



Non-Perspiring Knitting



Soft Knitting



Double Covered Thread



Non-Slip Band




The anti-bacterial properties of the lymphedema armsleeves produced by Aurafix's high quality double covered yarn and weaving technology have been developed and tested in accordance with the world standards. The yarns and knitting technique ensure maximum air permeability and minimizes perspiration by throwing out the moisture. Drip silicone elastic bands ensure that the product remains stable in the desired area and eliminates discomfort and uselessness caused by slip. They provide a friendly and non-irritating contact with the skin thanks to its soft knitting.





	Wrist to Shoulder Silicone Band	Wrist to Shoulder Shoulder Strap	Toe to Shoulder Silicone Band	Toe to Shoulder Shoulder Strap
CCL 1 18-21 mmHg	AV-130 ●	AV-150 ●	AV-133 ●	AV-153 ●
CCL 2 23-32 mmHg	AV-131 ●	AV-151 ●	AV-134 ●	AV-154 ●
CCL 3 34-46 mmHg	AV-132 ●	AV-152 ●	AV-135 ●	AV-155 ●

Size	cA	cC	cE	cG
1	15-16 cm	12-13 cm	17-22 cm	19-25 cm
2	16-19 cm	13-15 cm	20-26 cm	22-29 cm
3	18-21 cm	15-17 cm	23-29 cm	25-32 cm
4	20-23 cm	17-19 cm	26-32 cm	29-36 cm
5	22-25 cm	19-21 cm	29-35 cm	33-40 cm
6	24-27 cm	21-23 cm	32-38 cm	37-44 cm
7	26-29 cm	23-25 cm	35-41 cm	41-48 cm



Color



Caramel

YOUR SUPPORT IN LIFE





#yoursupportinlife

TRAVEL SOCKS

To stay still or sit in a confined space for long periods of time in long journeys over 4 hours, such as airplanes, cars, buses, or train journeys are unhealthy and may be dangerous.

Deep vein thrombosis is caused by blood obstruction and small clots. According to clinical studies conducted, travel socks reduce non-specific deep vein thrombosis and leg edema. Travel Socks apply a reduced compression towards the calf, with the highest pressure in the ankle. They reduce vessel diameter by gently compressing the veins and reduce blood clotting by increasing blood flow. Thus, travel socks prevent swelling in the legs, keep your legs healthy and energized, reduce the risk of deep vein thrombosis.

Aurafix Travel Socks, which have world-class pressure profile and knitting quality, are a comfortable solution with high quality and durable weaving. The raw materials used have the certificate of Oeko-Tex 100 and minimize allergy effects that may occur on the skin.



The anti-bacterial properties of the travel socks produced by Aurafix's high quality double covered yarn and weaving technology have been developed and tested in accordance with the world standards. The yarns and knitting technique ensure maximum air permeability and minimizes perspiration by throwing out the moisture. Thanks to its soft knitting, it provides a friendly and non-irritating contact with the skin. They ensure the product life to be longer by showing maximum resistance to pressure applied from the heel area.



Antibacterial Weaving



Non-Perspiring Knitting



Soft Knitting



Double Covered Thread



Pressure Resistant Heel



CCL 1
13-18 mmHg

A-D
Knee High

AV-125



Color




Black

Size

cB

- XS** 18-20 cm
- S** 20-23 cm
- M** 23-26 cm
- L** 26-29 cm
- XL** 29-32 cm



Ref VG-05



Compression Stockings Cladding Apparatus

Compression stockings are close knitted socks as they apply pressure on your veins in order to regulate your blood flow and can cause a compulsion especially while passing on your ankle.

Compression stockings cladding apparatus eliminates this condition and helps you wear your socks more quickly and easily. By means of this tool, you may prevent laddering, tearing or corrosion of your compression socks that have a delicate weaving and you prolong its life. They are used for all sizes including below knee, above knee (groin) and panty and all pressures. You can use them on all models in all compression socks with peep toe or closed.



SILVER
SOCKS



aurafix
SILVER SOCKS

GLOW WITH HEALTH
COMES FROM SILVER





I aurafix SILVER SOCKS

Aurafix Silver Socks, produced with pure silver technology, does not prevent blood circulation with an elastic band specially designed for people with diabetes, and regulates your blood circulation thanks to the silver knitting on the inner surface, thus preventing numbness and edema on your foot. With its antibacterial and anti-microbial properties, it prevents the wounds and infections that may occur in your feet, prevents the formation of bacteria that cause foot odor, and prevents the fungus and eczema diseases that may occur on your feet.

Thanks to the therapeutic properties of silver metal, they help to heal wounds and infections on your foot as soon as possible without festering. They isolate the static electricity accumulated in your body and cause you to get away from fatigue and stress. They prevent swelling of the foot which may occur due to various reasons and provide a feeling of comfort.



DIABETIC FEATURE

Aurafix Silver Socks have been specially designed for people with diabetes from pure silver yarn. Thanks to pure silver yarn technology, they keep your foot warm. Thanks to its antibacterial and antimicrobial properties, it prevents the wounds and infections that may occur in your feet. Due to the therapeutic properties of silver metal, they help to heal wounds and infections in your foot as soon as possible without festering.

Aurafix Silver Socks do not tighten your feet with the rubber specially designed for people with diabetes, they do not prevent blood circulation and they prevent the formation of edema thanks to the silver knitting on the inner surface. They do not irritate your feet because of the seamless nose feature, they prevent possible injuries.

FUNGUS AND ECZEMA PREVENTION FEATURE

They do not allow Silver bacteria to develop immunity. Bacteria need temperature, humidity and surfaces that prepare these conditions in order to survive. Therefore, socks create an ideal environment for bacteria. Silver removes these microorganisms and prevents them from forming colonies.

With antibacterial and antimicrobial effect that the pure silver yarn technology has, Aurafix Silver Socks prevent bacteria and infections that may occur on feet.

Thanks to their thermodynamic properties, they provide the sweat glands to function properly and provide a healthy usage by preventing the most important factor in the formation of sweat that is sweat.

It is recommended that Aurafix Silver Socks to be used pre-formed fungus and eczema as a treatment support product during the treatment process as well. This effect has been proven by clinical studies.





ANTI-ODOR FEATURE

Socks, produced with Aurafix pure silver yarn and natural fiber ensure that your foot let air in and thus prevent your feet from sweating. Silver yarn prevents the growth and reproduction of bacteria that spread sweat and foot odor.

Silver neutralizes sweat (ammonia) and deformed (denatured) proteins, eliminating odor molecules stemming from bacteria combination with its anti-bacterial properties.

ANTISTATIC FEATURE

Silver metal is the element having highest conductivity among known metals. Thanks to the anti-static feature of silver yarn technology, Aurafix Silver Socks reduce the feeling of fatigue by insulating the static electricity accumulated in the body and provide a more comfortable use.

They are also recommended as travel socks as they distribute the static electricity charge caused by air exchange during travel by means of silver ions.





THERMODYNAMIC FEATURE


Silver yarn reflects 95% of the heat energy coming out of your body back to its source and stores and distributes thermal energy. Thus, it keeps your feet warm in cold weather. Thanks to its thermal conductivity, the silver yarn distributes the heat in your feet to the outer environment in hot environments. So it keeps your feet cool. Silver yarn regulates your body's heat balance by ensuring that sweat glands work properly.

100% NATURAL

All of the products you see the Aurafix brand have a silver yarn made of pure silver. They do not contain any chemical and harmful substances. They are completely natural and anti-allergic. For this reason, all the properties added to the product with silver yarn technology do not lose its effect through the washing process. The specifications of the product are permanent throughout the life of the product.

They have been designed dedicatedly to ensure that the specifications of the product are permanent and the benefits of the silver are fully utilized.





Aurafix Silver Socks consist of a combination of pure silver, modal and polyamid. They do not contain synthetic yarn mixtures so that the maximum benefit from Silver Yarn technology is provided. Modal, produced from beech tree fibers, keeps moisture transfer and air permeability at the highest level while antimicrobial and antistatic properties of silver provide protection for your feet. The modal in the content not only supports abrasion resistance but also provides the texture to be soft and non-irritating. Thanks to the yarn quality used in our products, our products are friendly to your skin and long-lasting.



GLOW WITH HEALTH COMES FROM SILVER



SILVERSOCKS



#yoursupportinlife

Midcalf

DC-105



DC-105 Midcalf Silver Socks

Aurafix Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use them comfortably in your daily life with 3 color choices. They are suitable for unisex use.



Size	Number
S	35-38
M	39-42
L	43-46

Knee High

DC-120



DC-120 Knee High Silver Socks

Aurafix Knee High Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use below knee socks which are preferred especially during travelling and while doing sports comfortably in your daily life. Also suitable for unisex use.



Size	Number
S	35-38
M	39-42
L	43-46

Anklet

DC-110



DC-110 Anklet Silver Socks

Aurafix Anklet Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use ankle socks which are preferred especially while doing sports comfortably in your daily life with 4 color choices. They are suitable for unisex use.



Size	Number
S	35-38
M	39-42
L	43-46

Colors



.Black



Smoked



Navy



13%
Silver

Color



Black



25%
Silver

aurafix

Colors



Black



White



Smoked



Navy



13%
Silver

Midcalf Toe

DC-115



DC-115 Midcalf Toe Silver Socks

Aurafix Midcalf Toe Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. Especially in cases where your feet remain in the shoes too long as in the military service, they help keep your feet healthy by touching to your interdigital parts. They are suitable for unisex use.



Size	Number
S	35-40
M	41-44

Anklet Toe

DC-125



DC-125 Anklet Toe Silver Socks

Aurafix Anklet Toe Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use them comfortably in your daily life with black and white color choices. They are suitable for unisex use.



Size	Number
S	35-40
M	41-44

No-Show

DC-100



DC-100 No-Show Silver Socks

Aurafix No-Show Silver Socks are specially designed for patients with diabetes and regulate your blood circulation. In addition, they allow you to be protected from bacterial diseases such as foot odor and fungus. They do not tighten your foot with a specially manufactured elastic band; they do not prevent blood circulation. With the seamless feature, they do not irritate your feet. You can use them comfortably in your daily life with 4 color choices. They are suitable for unisex use.



Size	Number
S	35-38
M	39-42
L	43-46

Color



Black



13%
Silver

Colors



Black



White



10%
Silver

aurafix

Colors



Black



White



Smoked



Navy



12%
Silver

YOUR SUPPORT IN LIFE



PROMOTIONAL ITEMS



#yoursupportinlife

Posters




Compression Products 1

Ref PP-100



Compression Products 2

Ref PP-101



Lymphedema Armsleeves

Ref PP-102



Travel Socks

Ref PP-103



Silver Socks 1

Ref PP-800



Silver Socks 2

Ref PP-801

Floor Stickers



Compression Products 1

Ref PZ-100



Compression Products 2

Ref PZ-101



Silver Socks 1

Ref PZ-800



Silver Socks 2

Ref PZ-801

Roll Up I



Compression Products 1

Ref PR-100



Compression Products 2

Ref PR-101



Silver Socks

Ref PR-800

Standing Posters I



Compression Products

Ref PA-100



Silver Socks

Ref PA-800



Standing Mannequin

Ref PD-100



Compression Products Stand

Ref ST-02



Compression Products Prescription

Ref RECETE-VAR-EN



Carry Bags

Ref AC-S
45x36 cm

Ref AC-L
60x55 cm

Ref AC-M
45x36 cm



aurafix Point of Sale Sticker

Ref AS-EN



Below Knee Mannequins I



Above Knee & Arm Mannequins I



Full Body Mannequins



Follow Us



[/aurafix.tr](#)

aurafix
COMPRESSION PRODUCTS

YOUR SUPPORT IN LIFE

TSD Sağlık Ürünleri Tic. Ltd. Şti.
Adil Mah. Ekol Cad. No: 12/A 34935
Sultanbeyli - İSTANBUL/TURKEY
T:+90 216 399 30 55
F:+90 216 399 63 03
info@aurafix.com
aurafix.com

